

GetOn...with Life



DARING TO INNOVATE, FOR THE PHYSICALLY CHALLENGED.



ADAM MOTOR COMPANY

Adam Motor Company (AMC) is 100% wholly owned subsidiary of Omar Jibran Engineering Industries Ltd. (OJ).

Omar Jibran was incorporated in 1986. In Pakistan, Omar Jibran is one of the largest autoparts and other engineering products manufacturing company.

Adam Motor Company (AMC) was incorporated in 2001. Adam Motor Company (AMC) specializes in designing and manufacturing of mobility products. Realizing the vacuum in manufacturing of rehabilitation products for the ailing and the disabled, Adam Motor Company started a new division comprising the products like Wheelchairs both motorized and manual, Walking Sticks, Walkers, Stretchers and Hospital Beds.



QUALITY POLICY

"To manufacture and deliver each product and service to conform or exceed the customer specified quality on time at all times."



GetOn...*with Life*

Motorized Wheelchair




GetOn...with Life

Motorized Wheelchair

SALIENT FEATURES

- High quality steel.
- Dynamic electronic control.
- Stable front wheel drive.
- Independent drive wheel suspension
- Tyres Treads leaving no floor marks.
- For users weighting up to 150kg.
- High-back seat with headrest.
- Flip-up foot platform.
- Flexible location for mounting joystick.



SPECIFICATIONS:

Model No.	A1
Battery	12V x 2, 33AH
Motor Power	300W with Max. 200rpm
Controller	PG 50A
Electric Brake	Huafeng 24V
Overall Dimension	1130 x 640 x 1130mm
Weight w/o Battery	70kgs
Max. Speed (km/h)	7.5km/h
Turning Radius	0.55m
Max. Climbing Grade	12°
Wheel	260 x 85mm
Drive System	DG-198
Max. Loading Weight	100kgs
Total Weight w/bat.	77kgs
Ground Clearance	50 mm
Travel Range	up to 30km
Free Wheeling	Manual
Wheel Hub	Aluminum Alloy
Angle of Inclination:	10°
Max. Obstacle height:	90mm jump
Packing Size	90 x 69 x 79cm
Charger (included)	5A
Colors Available	Grey and Black

Warranty:

Frame 5 years, Battery, Electrical and Electronics 1 year.

GetOn...*with Life*

Foldable Motorized Wheelchair




GetOn...with Life

Foldable Motorized Wheelchair

SALIENT FEATURES

- High quality steel.
- Quickest folding & unfolding.
- Light weight body structure.
- Detachable Batteries.
- Carries upto 150kg or 330lbs.
- Universal fixed in-line joystick mount.
- Smooth & stable rare wheel drive.
- Five hours or 24km on each charge.
- Independent drive wheel suspension.
- Electronic controls.



SPECIFICATIONS:

Model No.	A2
Battery	24V 20AH
Motor Power	24V DC250W x 2
Controller	Intelligent Controller
Motor	2 x 250W Copper Brushed
Overall Dimension	1150 x 650 x 950mm
Weight w/o Battery	32kgs
Max. Speed (km/h)	0-6km/h (5 gears)
Turning Radius	0.55m
Max. Climbing Grade	≤ 13 °
Front Wheel	8" PU Tyre
Rear Wheel	12" Rubber Tyre
Max. Loading Weight	100kgs
Total Weight w/bat.	42kgs
Travel Range	10-20km
Base Width	45cm
Folding Size	800 x 400 x 660mm
Carton Size	850 x 430 x 760mm
Frame Material	High Carbon Steel
Angle of Inclination:	10°
Max. Obstacle height:	
Colors Available	Blue and Black

Warranty:

Frame 5 years, Battery, Electrical and Electronics 1 year.

GetOn...*with Life*

Lightweight Folding Manual Wheelchair



GetOn...with Life

Lightweight Folding Manual Wheelchair

SALIENT FEATURES

- High quality steel.
- Comfortable padded Armrest.
- Swing-out footrest to make transferring in & out easy.
- Carry pocket at the back.
- Convenient for Storage & transport with folding features.
- Provide comfort & mobility in an Economical Package.



SPECIFICATIONS:

Model No.	A3
Material	High Strength Steel
Front Wheel Size	8 inch
Rear Wheel Size	24 inch
Folding Size	900 x 900 x 280mm
Loading Capacity	140kg
Net Weight	14kg
Vehicle Width	680mm
Longest Height	900mm
Handrail Height	650mm
Handrail Length	350mm
Block Width	480mm
Seat Depth	450mm
Backrest Width	450mm

Warranty:
Frame one year.

GetOn...with Life


Walking Sticks



SALIENT FEATURES

- T Shaped Handle
- Z Shaped Handle
- Elbow Support Handle
- All walking sticks are made of lightweight anodized aluminum tube.
- The height of each stick is adjustable.
- The height can be adjusted by an easy to operate mechanism.
- The bottom of all sticks are provided with replaceable solid resin holder.

Various types of Walking Sticks to suit the need:



GetOn...with Life

Walkers



SALIENT FEATURES

- Designed to provide maximum support and confidence to the user.
- Easy-to-use, simple push-button mechanism for adjusting height and folding the walker.
- Each side folds independently to allow ease of movement through narrow spaces.

Various types of Walkers.



GetOn...with Life

RANGE OF REHABILITATION PRODUCTS

1. Motorized Wheelchair
2. Motorized Wheelchair (foldable)
3. Manual Wheelchairs
4. Walking Sticks & Crutches
5. Walkers
6. Patients Carrying Stretchers (foldable)
7. Hospital Beds Manual
8. Hospital Beds Motorized



Adam Motor Company Limited.

DSU-10, Pakistan Steel Downstream Industrial Estate Bin Qasim, Karachi, Pakistan.

Tel: +92 21 34723573-75 | Fax: +92 21 34723577

Email: sales@adammotor.pk | Website: adammotor.pk

